

# LifeView™ Patient Station

ati AmericanTeleCare



## Lead the Evolution. Connect Remotely with Patients using Video and Monitoring in one Device.

The American TeleCare® LifeView Video Patient Station combines the proven power of video telehealth and remote patient monitoring. This dynamic duo of patient monitoring and interactive video enables clinicians to monitor patient's health status, identify instability, and provide real time remote clinical intervention to re-stabilize patients. This evolutionary approach to health care finally delivers on the promise of providing the right care - at the right time - in the right place.

### The Power of Combined Video and Monitoring

- Achieve measurable outcomes as published
- Combine monitoring algorithms with face-to-face communications
- Visualize and assess patient status when their condition demands it
- Intervene with just-in-time education
- Utilize the movable camera to perform skin and wound assessments
- Address psycho-social needs
- Intervene at your patients most "vulnerable moments"

### Higher Quality Care at a Lower Cost

The LifeView Video Patient Station is unsurpassed in providing efficient healthcare in the home by reducing:

- Unscheduled physician office visits
- Skilled home care visits
- Hospitalizations/length of stay
- Emergency services

### Prepare for Health Care Evolution

The LifeView Patient Station offers healthcare organizations the ability to implement new strategic capabilities. Telehealth programs differentiate organizations while helping to improve patient compliance and increasing operational efficiency. Expand your business opportunities with:

- Disease management
- Private pay
- Hospice care giver support
- Palliative care
- Early discharge

### Created for the Patient

Designed with the patient in mind:

- **Ease of Readability:**
  - Large, high contrast, color screen
  - Large font, images and icons
- **Ease of Operation:**
  - Narrated instructions with answer verification
  - Multi-media instructions on use of peripheral medical devices
  - On-screen volume and navigation controls
- Compact size easily fits within patient's home
- Durable, to withstand life's little knocks
- Screen tilt adjustment reduces glare and enables patient use in sitting or standing position

## Powerful Features

Ideal for your most vulnerable patients, the LifeView Patient Station features two-way video/audio interface combined with medical peripherals and easy-to use monitoring algorithms. Consider the possibilities:

- Combines the features of the ATI inLife XP™ Patient Monitor with the audio-visual capabilities of the LifeView Patient Station
- Provides a critical view when needed for:
  - wound care
  - surgical incision
  - skin conditions
  - speech/physical therapy
  - medication verification
  - face to face teaching & support
  - psycho-social needs
  - edema
- Delivers the ability to listen to heart, lung and bowel sounds remotely
- Decreases false positives with vital signs taken in context
- Offers incomparable flexibility in scheduling and personalization of patient algorithms so you can deliver individualized care
- Delivers remote management of algorithm, vital sign, and software update capability - no need for a patient visit to make changes
- Maximizes ease of clinical management through the web based Provider Station software - results available to nurses and physicians from any computer browser, any where, any time
- Decreases patient teaching time due to superior design and ease-of-use
- Encrypted patient data is designed to operate within HIPAA regulations
- Provides text to speech capabilities

For detailed instructions in system operation, specifications, warnings, cautions, and additional information, please refer to the Operators Manual supplied with each product.

American TeleCare offers home telehealth programs tailored to your patient population and business models that include a variety of audio/video and monitoring stations with medical peripherals. Backed by superior customer service and technical support, American TeleCare solutions are at work helping to improve clinical outcomes, decrease costs, and increase clinician and patient satisfaction.

## We can show you real numbers, because we've seen real results

To arrange an assessment to explore the advantages of telehealth for your organization, **call 1-800-323-6667**

## Specifications

- Height: 10"; Depth: 9.7"; Weight: 26 lbs;
- Screen Size: 12.1"
- Connectivity: Analog POTS (Plain Old Telephone System) and broadband
- Power: Power supply input 100 -240 V AC, 50/60 Hz; 1.6A max 12 V DC
- Optional FDA Cleared Peripheral Medical Devices: Stethoscope, Scale, Blood Pressure Meter, Glucose Meter, Pulse Oximeter, Thermometer, PT/INR Meter
- USB and serial ports for connecting peripheral medical devices
- Designed with Microsoft Windows XP® Home operating system
- Integrated power surge protection
- 40 Gigabyte hard drive, 256K RAM, standard Intel 650 MHZ Celeron processor
- Movable camera

## Expansion Features

- HL-7 interface capable with other medical software
- Multi-media capability expansion for video-based training
- Custom connection to internet
- Multilingual text-to-speech
- CareTone® Telephonic Stethoscope - transfers clear, real-time heart and lung sounds to the Provider Station.



Practical solutions. Positive outcomes.

American TeleCare, Inc.  
15159 Technology Drive  
Eden Prairie, MN 55344-2273  
Customer and Clinical Services: 800-323-6667

Fax: 952-944-2247  
info@americantelecare.com  
Visit us at:

[www.americantelecare.com](http://www.americantelecare.com)